

# Purity Certificate

## Metagenics Purity Guarantee for Concentrated Omega-3 Formulas

<b>PRODUCT:</b>	CONCENTRATED OMEGA-3 ESSENTIAL FATTY ACIDS
<b>SOURCE:</b>	Selected, fresh, high EPA-DHA content saltwater fish
<b>DESCRIPTION:</b>	Third-party, purity-certified, expertly detoxified, high EPA-DHA marine lipid concentrate

Metagenics unconditionally guarantees the purity and quality of our concentrated omega-3 essential fatty acid formulas. For your peace of mind, we have crafted a line of omega-3 formulas that exceed the quality criteria established in the Council for Responsible Nutrition's Omega-3 Monograph.

Attention to every detail and exhaustive third-party and in-house testing enable us to make you the following promises:

### Metagenics concentrated omega-3 formulas are tested for:

Mercury	Less than 5 ppb (below detection limit)
Lead	Less than 25 ppb (below detection limit)
Dioxin	TEQ less than 5 ppt (below detection limit)
PCB's	TEQ less than 3 ppt (below detection limit)

- Also tested for yeasts, molds, and bacteria



### Metagenics concentrated omega-3 formulas contain:

- Only the highest quality, pharmaceutical grade omega-3 essential fatty acids
- EPA & DHA levels that meet label claim (milligram, *not just area percent*)
- Effective levels of antioxidant stabilizers

### Pharmaceutical grade

The manufacturing and purification technology used to produce the Metagenics line of concentrated omega-3 fatty acid formulas meets the exacting standards of the International Pharmacopoeia. These formulas are produced through an advanced molecular distillation process that ensures the highest purity and batch consistency.

### Gelatin softgels

The gelatin used in our softgels is obtained from countries known to be free of bovine spongiform encephalopathy (BSE) and have programs in place to prevent its occurrence.

### USP enteric coating

In order to avoid premature dissolution in the stomach, USP compliant enteric coating is utilized in some of our concentrated omega-3 formulas. USP guidelines stipulate that enteric coating must delay gastric disintegration by a minimum of one hour.

### Glossary of terms

**TEQ** – toxicity equivalency quotient; method used to represent the average toxicity potential when a type of toxin contains numerous subtypes (e.g., dioxins)

**PPB** – parts per billion

**PPT** – parts per trillion